

WIKI JS

Vickers, Knoop, Brinell




WIKI JS - THE MODELS



WIKI 100 JS

Motorized turret 6 positions for indenters and objectives.
 Motorized Z axis.
 Automatic measurement with autofocus.
 Software for automatic evaluation of indents.
 Automatic lighting.
 As optional: Manual XY stage 100x100 mm with 10 μ m step.

The automatic measurement minimizes subjectivity associated with human intervention. The tester can be used by different operators.



WIKI 200 JS

Top-of-the-line of automation and accuracy for Vickers and Knoop.
 Fully motorized system for Case Hardness Depth test **CHD**.
 Single or multiple samples in automatic cycle.
 Motorized turret 6 positions for indenters and objectives.
 Motorized Z axis.
 Automatic measurement with autofocus.
 Software for automatic evaluation of indents.
 Automatic lighting.
 Motorized XY stage 100x60 mm or 200x100 mm with 0,5 μ m step

Everything is automated, freeing users for other tasks and minimizing subjectivity associated with human intervention.

Standard Optional

WIKI JS LOAD FORCE RANGE

0.0098	0.0196	0.049	0.098	0.1471	0.1961	0.2452	0.4903	0.9807	1.961	2.942	4.903	9.807	19.61	29.42	49.03	98.07	196.1	294.2	490.3	980.7	N
0.001	0.002	0.005	0.01	0.015	0.02	0.025	0.05	0.1	0.2	0.3	0.5	1	2	3	5	10	20	30	50	100	kgf

WIKI JS 3 LOAD FORCE RANGE

0.0098	0.0196	0.049	0.098	0.1471	0.1961	0.2452	0.4903	0.9807	1.961	2.942	4.903	9.807	19.61	29.42	49.03	98.07	196.1	294.2	490.3	980.7	N
0.001	0.002	0.005	0.01	0.015	0.02	0.025	0.05	0.1	0.2	0.3	0.5	1	2	3	5	10	20	30	50	100	kgf

VICKERS FEASIBLE TESTS - DIN EN ISO 6507 / ASTM E-384

HV0.001	HV0.002	HV0.005	HV0.01	HV0.015	HV0.02	HV0.025	HV0.05	HV0.1	HV0.2	OHV.3	HV0.5	HV1	HV2	HV3	HV5	HV10	HV20	3HV0	HV50	HV100
---------	---------	---------	--------	---------	--------	---------	--------	-------	-------	-------	-------	-----	-----	-----	-----	------	------	------	------	-------

KNOOP FEASIBLE TESTS - DIN EN ISO 4545 / ASTM E-384

HK0.001	HK0.002	HK0.005	HK0.01	HK0.015	HK0.02	HK0.025	HK0.05	HK0.1	HK0.2
---------	---------	---------	--------	---------	--------	---------	--------	-------	-------

BRINELL DIN EN ISO 6506 / ASTM E-10 (Optional)

153.2	306.5	N	15.6	31.25	kgf	HBW2.5/15.6	HBW2.5/31.25
-------	-------	---	------	-------	-----	-------------	--------------

SUPERFICIAL R. DIN EN ISO 6508 / ASTM E-18 (Optional)

147.1	294.2	441.3	N	15	30	45	kgf
HR15 N/T/S/W/X/Y	HR30 N/T/S/W/X/Y	HR45 N/T/S/W/X/Y					

OPTIONAL TESTS:


Pack 1 HV/HK: HV0.001 - HV0.002 - HV0.005 (0.001 - 0.002 - 0.005 kgf)

Pack 2 HV: HV30 - HV50 (30 - 50 kgf)

Pack 3 HV: HV100 (100 kgf)

LOAD CELL AND CLOSED LOOP TECHNOLOGY

WIKI JS is the top-of-the-line of automation and accuracy for Vickers and Knoop hardness measurements in compliance with ASTM and ISO hardness standards. Load forces are applied through load cells and electronically controlled in "Closed Loop" (Pat. AFFRI) with a frequency of 1 khz. Each load force is automatically programmed and controlled assuring perfect linearity in every range eliminating the problems associated with traditional dead weight system testers. Results are not affected by any structural deflection, misalignment and external vibrations.



WIKI JS - MAIN FEATURES

MOTORIZED HEAD AND AUTOFOCUS




Up to 300 mm electronically controlled height capacity for fast or slow vertical movements. Very rapid and ultra-sensitive drive for a perfectly accurate autofocus. The autofocus combined with the automation of the whole software avoids human influence and gives repeatability even when used by different people.

6 SLOTS ROTATING TURRET



Horizontally rotating turret with four slots for magnification lenses and two for indenters. All optical microscope objectives can be pre-installed and combined with indenters for every Vickers and Knoop hardness scales. Optical objectives selection of 2.5x – 5x – 10x – 20x – 40x – 50x – 100x.

IMAGE AUTO-ANALYSIS



With software controlled focus, image cleaning, shading correction and regulated light source, reproducible results are obtained regardless of the number of indents measured. From perfectly polished to rough and etched samples, the auto-detection capabilities of WIKI JS allow measurements on a variety of sample surfaces.

LARGE AND STABLE BASE



The wide work table base is capable of bearing masses beyond 1000 kg which allows for steady hardness measurements on bulky or irregular pieces. It also offers a comfortable working base for small pieces.

X/Y MANUAL TABLE



Manual XY table 100x100 mm with 10 μ m step. This table is a perfect solution for not-daily multi indentation test cycles. The table allows manual CHD case depth tests and can be provided with digital micrometers for automatic CHD graph generation.

X/Y FULLY MOTORIZED



XY motorized table with an accuracy of ± 0.5 μ m steps. Reference points for indentation patterns can be positioned precisely where they are required. The table allows automatic multi-indentation CHD test cycles on multiple samples with perfect positioning on the entire area, no matter the indentations amount.

WIKI 100 JS SOFTWARE

High definition monitor 24"

Clean vision of the indent


Visual control of all results and live statistics

Direct conversion in HR, HB, HK and any other hardness scale.

Print results from template or save/import tests cycles from archive.


Control the whole instrument with the mouse.

Customizable test report with client logo, specimen information, statistics and graphs or export as CSV file.



The AFFRI Vickers measuring software has been studied to fulfill any client need and to be accessible to every operator. This is a **"SMART SOFTWARE"** which results extremely easy to be used and can be customized to display only needed testing procedures.

ONLY 5 ICONS TO GET RESULTS:



1 - MAGNIFICATION



2 - AUTOFOCUS



3 - AUTO-LIGHT




4 - TEST METHOD




5 - AUTOMATIC MEASURE

SMART SOFTWARE - MAIN FEATURES


REAL MAGNIFICATION - Thanks to the motorized turret, different lenses can be selected with a simple click. Digital zoom is also available.




AUTOFOCUS - Motorized focus is always at the right linear quote.




AUTO-LIGHTING - Automatic light regulation on any surface.




TEST METHOD SELECTION - Only one window for the selection of everything you need for the test.




AUTO MEASURE ON CRITICAL SURFACES - From perfectly polished to rough & etched samples, the software will automatically measure indents on any sample surface.



MANUAL MEASURE - Manual indent evaluation.



AUTOMATIC CONVERSION - From standards to hardness scale.




DYNAMIC RESULTS

Color highlighted results and live statistics. Watch result list and edit or modify single tests.



LIVE GRAPHS

Choose between 4 graphs. Print results from template or save and import tests cycles from archive.



WIKI 200 JS SOFTWARE




Get real and detailed images of each sample.

Visual control of each sample and patterns.

See a list of each sample and each pattern.
Save or import traverses, edit, move, copy and paste.

Control the whole instrument with the mouse including dynamic movements on X/Y/Z axis at fast or slow speed.



Customizable test report with client logo and information. Put as many sample information as you need, including pictures of the specimen and test area. See the pattern and each indentation at a glance. Examine results, statistics and CHD diagram with outlined depth. Results can be exported also as CSV file.

SINGLE OR MULTIPLE SAMPLES IN AUTOMATIC CYCLE: Just map out indentation traverses where they are required, set the load and press START, the hardness tester intelligently follows the predefined patterns, indents the sample, focuses when needed, measures, and generates data dynamically.

The software is designed for an intuitive and simple use. With three easy steps, it provides added precision when positioning indents thanks to its integrated macro view technique and layout tools. By visualizing the complete sample or a single sample, traverses and/or patterns can now be mapped out with unequalled precision.

01 SURFACE MAPPING

Obtain a perfect, detailed, high resolution view of the whole sample holder offering sharp close-ups as well as global views. Multi sample vision offers a complete image of a sample without any distortion, no matter its size. No need for a second camera.



02 SETUP PATTERNS

Save, open, modify, copy and paste or create new patterns to predefined locations with a simple click. Traverses and patterns can be individually adjusted. Create case depth traverses or fill a surface with indentation points to control sample uniformity.



Samples can be mapped separately. Singular identification allows to group patterns regarding only one sample. Results storing and reporting will be well organized.


Traverse's indents can be modified at any time: delete, move or place a new one. Set where to perform autofocus and choose a different test load for each pattern if needed.

Rotate traverses with T-Bar tool.

Reports are automatically generated and archived during the test cycle. Just click on one plot to check the indentation.

03 PRESS START GET RESULTS

The software follows the patterns, indents the sample, measures, and generates data dynamically. Review results in graphical and/or tabular format. Export results to any spreadsheet application, or simply print standard or customized reports.



EVERYTHING IS AUTOMATED, FREEING USERS FOR OTHER TASKS: Auto focusing, automatic measuring and reporting, allows this system to function unattended for hours without interruption, saving time and money, thus increasing output and productivity.

AUTOMATIC MEASUREMENT CYCLE ON MULTIPLE SAMPLES

CHD TRAVERSES AND PATTERNS

Single or multiple traverses/patterns can be rapidly created. With one simple click of the mouse the line, the angle and the starting point of the indentations are determined.

CREATE, SAVE AND RELOAD

No need to create the same pattern over and over again. This feature is extremely useful for users who analyze the same kind of areas repeatedly. Once a pattern has been created, you can save it and re-load it later to duplicate the analysis on a new sample.

EDIT, MOVE OR DELETE


Select the pattern from preview image and modify direction, position, hardness scale, number and distances of indent, pattern name and point of focus.

ZOOM IN AND OUT

Zoom out to look at the entire sample in order to identify the pattern position and direction. Zoom in to verify pattern distances and spot surface imperfections.

Traverse layouts and test points can be programmed by simply clicking on the desired test point locations.


Save, cut, copy, or paste traverses/patterns to predefined locations with a simple click of the mouse.




Use image preview to control the whole sample patterns. Identify the test zone and place patterns with a simple mouse click.

The T-Bar tool rotates traverses to any angle to ensure its perpendicularity to the sample edge or to accommodate sample tilts.


CHD TRAVERSES
Case Hardness Depth




FREE PATTERNS
Create the shape you want



SURFACE PATTERN
Study sample uniformity




RADIUS PATTERN
Study uniformity per radius



LENGTH AND DRAWING TOOLS

Use this feature to place reference points, find the middle, text some notes and measure the length of anything.



No matter the complexity of the pattern layout, combining intelligent macro views with drawing and length tools allows traverses and/or patterns to be positioned precisely where they are needed.

These features also help finding the center of the sample or the center of a gear tooth flank and the pitch point. Placing traverses and pattern is made easy.

Use this tool to draw straight and/or parallel lines, add text notes and measure lengths for report purposes.

RESULTS INSTANTANEOUS DATA REVIEW

Following an automated run, individual indents can be tracked by clicking on the numbered impression. Intelligent software accurately remembers where the impression was made and automatically moves the stage to the chosen indent.

You can choose to not include, re-measure the impression manually with the movable gridlines or make a new indent. When excluded or re-measured, statistics are updated on the fly. Instant graphical view of Effective Case Depth.

REPORT CREATION

Print results directly from the software or export data to the spreadsheet program of your choice for further statistical analysis.


Select the sample and the pattern.

Spot tested zone from number labels.

See each result including depth, diagonals and HRC conversion or others.

See dynamic statistics, define CHD, RHT, NHT and three custom case depths.

See Case Depth chart where hardness value is plotted against depth. Identify specified case depth values from dotted lines.



Export single or all reports in PDF.

Export data as text file or CSV.

See preselected statistics including max, min, mean and deviation.

See uniformity chart with high and low highlighted tolerances.

Detect unexpected results, click the plot to go to indent. Verify, measure again or replace with a new indent.


AUTOMATIC MEASUREMENT CYCLE ON WELDS

TESTING ON WELDS


After scanning the whole sample, the fusion zones and the HAZ are clearly visible and distinguishable from the base material, even with a scratched surface.

Using software's tools it is easy to draw different indentation patterns lines with correct and precise positioning. In less than one minute add pattern positions with determined spacing between indents, defined distance from the border, from inside and outside surface, from fusion line or weld centerline.

The entire indentation distance can be simply measure in one complete view and added to your final report.





Scanned images of weld cross section.



Customizable test report with client logo and sample information. See weld image, pattern, results and statistics at a glance.

See hardness curves grouped all together.



REFERENCE CIRCLE TOOL

Ideal for irregular or curved samples, where indents need to be at a given distance from the edge. This exceptional tool allows indents to be positioned at precise distances from the sample's edge. Once a radius is specified, the software shows a red circle around the indent position. Use this reference to go across a border and/or zone and add the indent position.

This visual guide is the fastest way to create a defined pattern on a multi-zoned or irregular sample.


CUSTOM STEP TOOL

When the tool is active, the virtual Joystick commands move the stage for the distance defined by step instead of a continuous movement.

Use this feature to create a pattern line where groups of spaced indents are placed in base material, HAZ and fusion. Define and fix a custom stage movement in X, Y or Z axis. Move on the point of interest (e.g. fusion line) and add an indent position for each spacing distance step.

REAL TIME SUPPORT:

The remote control connects AFFRI's testers from anywhere in the world with AFFRI's engineers. Our experts can remotely diagnose any technical issue, provide additional operator training and update the software version.



CUSTOMIZABLE REPORT TEMPLATES

Report are created directly by Affri team. Templates can be 100% customized basing on any need. Not only company logo, sample map and indentation images, the report layout can be filled with any test information. Time by time, after the test cycle, choose what to export in final report by flagging or unflagging simple options.

SL	NAME	TEST NAME	TEST NUMBER	VALUE	UNIT	SL	NAME	TEST NAME	TEST NUMBER	VALUE	UNIT
1	Sample 1	Hardness	1	100	MPa	1	Sample 1	Hardness	1	100	MPa
2	Sample 2	Hardness	2	120	MPa	2	Sample 2	Hardness	2	120	MPa
3	Sample 3	Hardness	3	150	MPa	3	Sample 3	Hardness	3	150	MPa
4	Sample 4	Hardness	4	180	MPa	4	Sample 4	Hardness	4	180	MPa
5	Sample 5	Hardness	5	200	MPa	5	Sample 5	Hardness	5	200	MPa
6	Sample 6	Hardness	6	220	MPa	6	Sample 6	Hardness	6	220	MPa
7	Sample 7	Hardness	7	250	MPa	7	Sample 7	Hardness	7	250	MPa
8	Sample 8	Hardness	8	280	MPa	8	Sample 8	Hardness	8	280	MPa
9	Sample 9	Hardness	9	300	MPa	9	Sample 9	Hardness	9	300	MPa
10	Sample 10	Hardness	10	320	MPa	10	Sample 10	Hardness	10	320	MPa
11	Sample 11	Hardness	11	350	MPa	11	Sample 11	Hardness	11	350	MPa
12	Sample 12	Hardness	12	380	MPa	12	Sample 12	Hardness	12	380	MPa
13	Sample 13	Hardness	13	400	MPa	13	Sample 13	Hardness	13	400	MPa
14	Sample 14	Hardness	14	420	MPa	14	Sample 14	Hardness	14	420	MPa
15	Sample 15	Hardness	15	450	MPa	15	Sample 15	Hardness	15	450	MPa
16	Sample 16	Hardness	16	480	MPa	16	Sample 16	Hardness	16	480	MPa
17	Sample 17	Hardness	17	500	MPa	17	Sample 17	Hardness	17	500	MPa
18	Sample 18	Hardness	18	520	MPa	18	Sample 18	Hardness	18	520	MPa
19	Sample 19	Hardness	19	550	MPa	19	Sample 19	Hardness	19	550	MPa
20	Sample 20	Hardness	20	580	MPa	20	Sample 20	Hardness	20	580	MPa
21	Sample 21	Hardness	21	600	MPa	21	Sample 21	Hardness	21	600	MPa
22	Sample 22	Hardness	22	620	MPa	22	Sample 22	Hardness	22	620	MPa
23	Sample 23	Hardness	23	650	MPa	23	Sample 23	Hardness	23	650	MPa
24	Sample 24	Hardness	24	680	MPa	24	Sample 24	Hardness	24	680	MPa
25	Sample 25	Hardness	25	700	MPa	25	Sample 25	Hardness	25	700	MPa
26	Sample 26	Hardness	26	720	MPa	26	Sample 26	Hardness	26	720	MPa
27	Sample 27	Hardness	27	750	MPa	27	Sample 27	Hardness	27	750	MPa
28	Sample 28	Hardness	28	780	MPa	28	Sample 28	Hardness	28	780	MPa
29	Sample 29	Hardness	29	800	MPa	29	Sample 29	Hardness	29	800	MPa
30	Sample 30	Hardness	30	820	MPa	30	Sample 30	Hardness	30	820	MPa
31	Sample 31	Hardness	31	850	MPa	31	Sample 31	Hardness	31	850	MPa
32	Sample 32	Hardness	32	880	MPa	32	Sample 32	Hardness	32	880	MPa
33	Sample 33	Hardness	33	900	MPa	33	Sample 33	Hardness	33	900	MPa
34	Sample 34	Hardness	34	920	MPa	34	Sample 34	Hardness	34	920	MPa
35	Sample 35	Hardness	35	950	MPa	35	Sample 35	Hardness	35	950	MPa
36	Sample 36	Hardness	36	980	MPa	36	Sample 36	Hardness	36	980	MPa
37	Sample 37	Hardness	37	1000	MPa	37	Sample 37	Hardness	37	1000	MPa
38	Sample 38	Hardness	38	1020	MPa	38	Sample 38	Hardness	38	1020	MPa
39	Sample 39	Hardness	39	1050	MPa	39	Sample 39	Hardness	39	1050	MPa
40	Sample 40	Hardness	40	1080	MPa	40	Sample 40	Hardness	40	1080	MPa
41	Sample 41	Hardness	41	1100	MPa	41	Sample 41	Hardness	41	1100	MPa
42	Sample 42	Hardness	42	1120	MPa	42	Sample 42	Hardness	42	1120	MPa
43	Sample 43	Hardness	43	1150	MPa	43	Sample 43	Hardness	43	1150	MPa
44	Sample 44	Hardness	44	1180	MPa	44	Sample 44	Hardness	44	1180	MPa
45	Sample 45	Hardness	45	1200	MPa	45	Sample 45	Hardness	45	1200	MPa
46	Sample 46	Hardness	46	1220	MPa	46	Sample 46	Hardness	46	1220	MPa
47	Sample 47	Hardness	47	1250	MPa	47	Sample 47	Hardness	47	1250	MPa
48	Sample 48	Hardness	48	1280	MPa	48	Sample 48	Hardness	48	1280	MPa
49	Sample 49	Hardness	49	1300	MPa	49	Sample 49	Hardness	49	1300	MPa
50	Sample 50	Hardness	50	1320	MPa	50	Sample 50	Hardness	50	1320	MPa
51	Sample 51	Hardness	51	1350	MPa	51	Sample 51	Hardness	51	1350	MPa
52	Sample 52	Hardness	52	1380	MPa	52	Sample 52	Hardness	52	1380	MPa
53	Sample 53	Hardness	53	1400	MPa	53	Sample 53	Hardness	53	1400	MPa
54	Sample 54	Hardness	54	1420	MPa	54	Sample 54	Hardness	54	1420	MPa
55	Sample 55	Hardness	55	1450	MPa	55	Sample 55	Hardness	55	1450	MPa
56	Sample 56	Hardness	56	1480	MPa	56	Sample 56	Hardness	56	1480	MPa
57	Sample 57	Hardness	57	1500	MPa	57	Sample 57	Hardness	57	1500	MPa
58	Sample 58	Hardness	58	1520	MPa	58	Sample 58	Hardness	58	1520	MPa
59	Sample 59	Hardness	59	1550	MPa	59	Sample 59	Hardness	59	1550	MPa
60	Sample 60	Hardness	60	1580	MPa	60	Sample 60	Hardness	60	1580	MPa
61	Sample 61	Hardness	61	1600	MPa	61	Sample 61	Hardness	61	1600	MPa
62	Sample 62	Hardness	62	1620	MPa	62	Sample 62	Hardness	62	1620	MPa
63	Sample 63	Hardness	63	1650	MPa	63	Sample 63	Hardness	63	1650	MPa
64	Sample 64	Hardness	64	1680	MPa	64	Sample 64	Hardness	64	1680	MPa
65	Sample 65	Hardness	65	1700	MPa	65	Sample 65	Hardness	65	1700	MPa
66	Sample 66	Hardness	66	1720	MPa	66	Sample 66	Hardness	66	1720	MPa
67	Sample 67	Hardness	67	1750	MPa	67	Sample 67	Hardness	67	1750	MPa
68	Sample 68	Hardness	68	1780	MPa	68	Sample 68	Hardness	68	1780	MPa
69	Sample 69	Hardness	69	1800	MPa	69	Sample 69	Hardness	69	1800	MPa
70	Sample 70	Hardness	70	1820	MPa	70	Sample 70	Hardness	70	1820	MPa
71	Sample 71	Hardness	71	1850	MPa	71	Sample 71	Hardness	71	1850	MPa
72	Sample 72	Hardness	72	1880	MPa	72	Sample 72	Hardness	72	1880	MPa
73	Sample 73	Hardness	73	1900	MPa	73	Sample 73	Hardness	73	1900	MPa
74	Sample 74	Hardness	74	1920	MPa	74	Sample 74	Hardness	74	1920	MPa
75	Sample 75	Hardness	75	1950	MPa	75	Sample 75	Hardness	75	1950	MPa
76	Sample 76	Hardness	76	1980	MPa	76	Sample 76	Hardness	76	1980	MPa
77	Sample 77	Hardness	77	2000	MPa	77	Sample 77	Hardness	77	2000	MPa
78	Sample 78	Hardness	78	2020	MPa	78	Sample 78	Hardness	78	2020	MPa
79	Sample 79	Hardness	79	2050	MPa	79	Sample 79	Hardness	79	2050	MPa
80	Sample 80	Hardness	80	2080	MPa	80	Sample 80	Hardness	80	2080	MPa
81	Sample 81	Hardness	81	2100	MPa	81	Sample 81	Hardness	81	2100	MPa
82	Sample 82	Hardness	82	2120	MPa	82	Sample 82	Hardness	82	2120	MPa
83	Sample 83	Hardness	83	2150	MPa	83	Sample 83	Hardness	83	2150	MPa
84	Sample 84	Hardness	84	2180	MPa	84	Sample 84	Hardness	84	2180	MPa
85	Sample 85	Hardness	85	2200	MPa	85	Sample 85	Hardness	85	2200	MPa
86	Sample 86	Hardness	86	2220	MPa	86	Sample 86	Hardness	86	2220	MPa
87	Sample 87	Hardness	87	2250	MPa	87	Sample 87	Hardness	87	2250	MPa
88	Sample 88	Hardness	88	2280	MPa	88	Sample 88	Hardness	88	2280	MPa
89	Sample 89	Hardness	89	2300	MPa	89	Sample 89	Hardness	89	2300	MPa
90	Sample 90	Hardness	90	2320	MPa	90	Sample 90	Hardness	90	2320	MPa
91	Sample 91	Hardness	91	2350	MPa	91	Sample 91	Hardness	91	2350	MPa
92	Sample 92	Hardness	92	2380	MPa	92	Sample 92	Hardness	92	2380	MPa
93	Sample 93	Hardness	93	2400	MPa	93	Sample 93	Hardness	93	2400	MPa
94	Sample 94	Hardness	94	2420	MPa	94	Sample 94	Hardness	94	2420	MPa
95	Sample 95	Hardness	95	2450	MPa	95	Sample 95	Hardness	95	2450	MPa
96	Sample 96	Hardness	96	2480	MPa	96	Sample 96	Hardness	96	2480	MPa
97	Sample 97	Hardness	97	2500	MPa	97	Sample 97	Hardness	97	2500	MPa
98	Sample 98	Hardness	98	2520	MPa	98	Sample 98	Hardness	98	2520	MPa
99	Sample 99	Hardness	99	2550	MPa	99	Sample 99	Hardness	99	2550	MPa
100	Sample 100	Hardness	100	2580	MPa	100	Sample 100	Hardness	100	2580	MPa

TABULAR AND TXT DATA EXPORT

Print the results directly from the Affri software or export data in txt, csv and tabular format to the spreadsheet program of your choice for further statistical analysis. Images and statistics can be saved or copied easily and laid out in a standard, or customized, MS Office templates. Use "paste as link" native tool to update your data sheet with exported results and build a well organized database.

ACCESSORIES

Affri provides a large variety of accessories to fulfill any purpose of test. Customized solution based on your needs can be made for perfect tests on rough pieces. A series of different anvils is available to test every size of test piece. Variety of accessories to facilitate testing on small or oddly shaped items. All AFFRI's accessories are customizable according to customers specifications, depending on dimensions and geometry of the samples and finished products.



TEST BLOCKS

Micro Vickers Art. A004.0.008
Knoop Art. A004.0.010
Test blocks with specific values are available.



OBJECTIVES

2,5X W001.0.006 / 5X W001.0.000 /
10X W001.0.001 / 20X W001.0.002 /
40X W001.0.003 / 50X W001.0.004 /
100X W001.0.005




BENCH SUPPORT TABLE

Side table suitable right or left
Art. A010.0.024



INDENTERS

Vickers Art. 700.1.5.027
Knoop Art. 700.1.5.028
Dual indenter shaft Art. 700.1.5.029



MOTORISED TABLE 200x200mm

Travel 200x100mm with 0,5µm step
Art. A055.0.001



MOTORISED TABLE 150x150mm

Travel 100x60mm with 0,5µm step
Art. A055.0.002



MANUAL TABLE

100x100mm with 10µm step
Art. A009.0.001



CLAMPING VICE

Adjustable from 0 to 50mm
Art. A049.1.001



SINGLE SAMPLE HOLDER

Self level sample holder
Art. A.055.0.006 (Insert ring is needed)
Art. A.055.0.014 (All sample diameter)



MUTLI SAMPLE HOLDER

Up to 10 samples per time
Art. A055.0.003 (Base 200x100 mm)
Art. A055.0.004 (Drawer 10x30mm)
Art. A055.0.005 (Drawer 8x40mm)
Other drawer are available



EXCLUSIVE DESIGN

The innovative design of WIKI JS is AFFRI's unique and exclusive. Comfortable and ergonomic working station built for facilitate operator's movements, allowing an organized and well-ordered work. Built-in side case to preserve tester's accessories.



WORKING STATION

Solid and compact workbench with a large locker to accommodate computer, console and cables. Wide side table, suitable right or left, to have all within operator's reach.

WIKI JS

FORCE RANGE

Vickers/Knoop:	0.0098 - 0.0196 - 0.049 - 0.098 - 0.1471 - 0.1961 - 0.2452 - 0.4903 - 0.9807 - 1.961 - 2.942 - 4.903 - 9.807 - 19.61 - 29.42 - 49.03 - 98.07 - 196.1 - 294.2 - 490.3 - 980.7 N (0.001 - 0.002 - 0.005 - 0.01 - 0.015 - 0.02 - 0.025 - 0.05 - 0.1 - 0.2 - 0.3 - 0.5 - 1 - 2 - 3 - 5 - 10 - 20 - 30 - 50 - 100 kgf)
Brinell:	153.2 - 306.5 N (15.6 - 31.25 kgf)
Superficial Rockwell:	147.1 - 294.2 - 441.3N (15 - 30 - 45 kgf)

WIKI 100/200 JS FEASIBLE TESTS

Vickers:	HV0.01 - HV0.015 - HV0.02 - HV0.025 - HV0.05 - HV0.1 - HV0.2 - HV0.3 - HV0.5 - HV1 - HV2 - HV3 - HV5 - HV10
Knoop:	HK0.01 - HK0.015 - HK0.02 - HK0.025 - HK0.05 - HK0.1 - HK0.2 - HK0.3 - HK0.5 - HK1 - HK2

WIKI 100/200 JS 3 FEASIBLE TESTS

Vickers:	HV0.1 - HV0.2 - HV0.3 - HV0.5 - HV1 - HV2 - HV3 - HV5 - HV10 - HV20 - HV30
Knoop:	HK0.1 - HK0.2 - HK0.3 - HK0.5 - HK1 - HK2

OPTIONAL TESTS (Applicable to all models)

Superficial Rockwell:	HR15N - HR30N - HR45N - HR15T - HR30T - HR45T - HR15S - HR30S - HR45S - HR15W - HR30W - HR45W - HR15X - HR30X - HR45X - HR15Y - HR30Y - HR45Y
Brinell HBW / HBWT:	2.5/15,6 - 2.5/31.25

TECHNICAL DATA

Accuracy:	Better than 0.1 %
Principle of Operation:	Load Cell and Closed Loop (Affri patent)
Standards:	EN-ISO 6506 / EN-ISO 6507 / EN-ISO 6508 / ASTM-E384 / EN-ISO 4545 / ASTM-E92 / ASTM E10 / ASTM E08 / ASTM E103 / JIS
Vertical Stroke:	Motorized 240 mm / 9.4" (as optional 300 mm / 12" or 700 mm / 27.5")
Depth Capacity:	135 mm / 5.5"
Turret:	Automatic and motorized - 6 positions
Indenter:	Vickers - As option Knoop and Brinell
Camera:	1.3 MP USB2 B/W HD
Focus and Reading:	Automatic and manual
Lighting:	Energy Efficient Cool LED Light Source
Network:	Wire connection for technical assistance and auto-diagnosis
X-Y Table:	Motorized with 0.5 ?m steps 100 x 60 mm / 3.9 x 2.3" or 200 x 100 mm / 7.8 x 3.9"
Tolerable weight:	50 kg
Dwell Time:	From 5 to 60 seconds programmable
Temperature Range:	From 10 °C to 35 °C
Data Output:	USB / Ethernet
Power Supply:	110 or 220 V / 50÷60 Hz
Software:	Affri - OMAG
Fields Of Use:	For micro and macro Vickers and case depth test on every metals: iron, steel, tempered steel, cast iron, brass, aluminium, copper and metal alloys. Heat treatment, hardening, nitriding, cementation and hardfacing. Knoop test on ceramic and glass materials.
Packaging:	120 x 120 x 160 cm / 47 x 47 x 65" - 160/200 kg



Made by:
OMAG di AFFRI D. S.r.l.
 Via M. Tagliaferro, 8, I-21056 INDUNO OLONA - CEE (VA) - ITALY
 Tel. +39 0332 200546 Fax +39 0332 203704
 info@omagaffri.com

Europe/Asia:
AFFRI®
 Via M. Tagliaferro, 8, I-21056 INDUNO OLONA - CEE - (VA) - ITALY
 Tel. +39 0332 201533 +39 0332 206289 Fax +39 0332 203621
 info@affri.com - www.affri.com

America:
AFFRI Inc.
 850 Dillon Dr. Wood Dale, IL 60191
 Tel. 224 374 0931 - 630 303 1588
 sales@affriusa.com - www.affri.com